# OHAUS Label Printing Instruction Manual

English

Español

Français

Deutsch

Italiano

Polski

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Operation Videos

EX Label Printing Function

OHAUS Label Designer

## Table of Contents

<b>1. INTRODUCTION</b> .....	<b>2</b>
<b>2. LABEL PRINTING FUNCTION</b> .....	<b>2</b>
2.1 Explorer .....	2
2.1.1 Step 1: Connect Explorer to label printer .....	2
2.1.2 Step 2: Enable label printing function in Explorer .....	3
2.1.3 Step 3: Print labels .....	3
2.2 Pioneer PX .....	4
2.2.1 Step 1: Connect Pioneer PX to label printer .....	4
2.2.2 Step 2: Enable label printing function in Pioneer PX .....	5
2.2.3 Step 3: Print labels .....	5
<b>3. PC SOFTWARE - OHAUS LABEL DESIGNER</b> .....	<b>6</b>
3.1 Where to download the OHAUS Label Designer .....	6
3.2 How to use OHAUS Label Designer .....	6
3.2.1 Step 1: Connect balance to PC .....	6
3.2.2 Step 2: Startup .....	7
3.2.3 Step 3: Customize Template .....	9
3.2.4 Step 4: Write label template to balance .....	13

## 1. INTRODUCTION

This manual is to help set up and configure an OHAUS Explorer balance (software version 2.23 and later) or a Pioneer PX balance (software version 1.08 and later) for printing labels and weighing results by a ZPL-compatible label printer. Please read this manual completely before software installation.

This label printing function of Explorer/Pioneer balance has been tested with a Zebra GC420t. It can be compatible with any label printer, which uses the ZPLII programming language.

As label formats can be highly customized, this manual is not intended as a tutorial on how to customize a label printed by the label printer.

## 2. LABEL PRINTING FUNCTION

### 2.1 Explorer


Explorer series balance allows print labels in 3 steps.


- **Step 1:** Connect the balance to a label printer using an RS232 cable
- **Step 2:** Enable label printing function in Explorer
- **Step 3:** Print labels by pressing **Print** in Explorer

#### 2.1.1 Step 1: Connect Explorer to label printer

For connecting an OHAUS Explorer balance to a label printer, you need to prepare the following items:

- A male-to-male RS232 cable (OHAUS PN 30057595)
- An OHAUS Explorer balance (software version 2.23 and later)
- A label printer that supports ZPL commands

**Note:**

Usually, Explorer can be connected to a label printer with default communication settings as right.

Before printing, make sure that both balance and printer have the same settings i.e., Baud Rate, Transmission, and Handshake.

To change the connection settings of balance, enter “Menu – Communication – RS-232 Standard – Device Settings”.

## 2.1.2 Step 2: Enable label printing function in Explorer

**Label Printing** menu is available in Explorer balance (software version 2.23 and later) with default **Off**. To enable label printing, enter “Menu – Communication – RS-232 Standard” and set the **Label Printing** menu to **ON**.

Then **Label Template** menu appears with default English label template “100mm \* 75mm”.

## 2.1.3 Step 3: Print labels

To print labels, choose label template and press **Print** button in Explorer.

There are 2 built-in templates in Explorer.

- #1 - 100mm \* 75mm
- #2 - 100mm \* 10mm

Below are **printout examples** using default templates:

Template	#1 - 100mm * 75mm	#2 - 100mm * 10mm
Size		
Content	<ul style="list-style-type: none"> <li>• Data/Time</li> <li>• Sample ID</li> <li>• Sample Name</li> <li>• Batch ID</li> <li>• Result</li> <li>• Gross Weight</li> <li>• Tare Weight</li> <li>• Net Weight</li> <li>• User Name</li> </ul>	<ul style="list-style-type: none"> <li>• Result</li> </ul>

When Explorer is in label printing mode, it prints stable results only.

Please refer to **Section 3 PC Software - OHAUS Label Designer** for instruction to change the default label templates.

## 2.2 Pioneer PX

### 2.2.1 Step 1: Connect Pioneer PX to label printer

Use this RS232 interface to connect the balance to a printer	Male-to-male RS232 cable

**Note:**

Usually, Pioneer PX can be connected to a label printer with default communication settings (baud rate **9600**, Transmission **8-No-1**, Handshake **None**).

Before printing, make sure that both balance and printer have the same settings i.e., Baud Rate, Transmission, and Handshake.

To change the connection settings of balance, long press **Menu** and enter **RS232** menu.

**2.2.2 Step 2: Enable label printing function in Pioneer PX**

**Label Printing** menu is available in Pioneer PX balance (software version 1.08 and later) with default **Off**. To enable label printing, enter “Menu – Print – Label Printing” and set the **Label Printing** menu to **ON**.

Pioneer PX balance comes with one default English label template.

**2.2.3 Step 3: Print labels**

To print labels, press **Print** button in Pioneer PX.

There are only one built-in template. Below is the **printout example** using default template:

<b>Size</b>	<p style="text-align: center;">← 100 mm →</p> <p style="text-align: right;">75 mm</p>
<b>Content</b>	<ul style="list-style-type: none"> <li>• Data/Time</li> <li>• Balance Name</li> <li>• Balance ID</li> <li>• Gross Weight</li> <li>• Tare Weight</li> <li>• Net Weight</li> <li>• Signature</li> </ul>

When Pioneer PX is in label printing mode, it prints stable results only.

Please refer to **Section 3 PC Software - OHAUS Label Designer** for instruction to change the default label templates.

### 3. PC SOFTWARE - OHAUS LABEL DESIGNER

OHAUS Label Designer enables the user to edit label templates and save to Explorer/Pioneer PX by connecting to a computer.

#### 3.1 Where to download the OHAUS Label Designer

Click to get free [OHAUS Label Designer](#) software.

You can also download the latest software from the official website **OHAUS.com** at the path:

#### 3.2 How to use OHAUS Label Designer

##### 3.2.1 Step 1: Connect balance to PC

This software works with the software version of Explorer balance (software version 2.23 and later) or Pioneer PX balance (software version 1.08 and later).

The operating system of PC for running OHAUS Label Designer is Windows 10® 32 / 64 bit.

To connect the balance with PC, you need the below accessory.

Interface kit RS232-USB (OHAUS PN 30304101)



### 3.2.2 Step 2: Startup




OHAUS Label Designer allows connecting balance to PC automatically or manually.

When connecting to Explorer, ensure the balance has been logged in as Admin or Supervisor before running the OHAUS Label Designer software.

#### Automatic Connection

 <p style="text-align: center;">Figure 3.2.2-1</p>	<ol style="list-style-type: none"> <li>1. Double-click the software.</li> <li>2. Choose the operation Language in the upper right corner. Below languages are available. <ul style="list-style-type: none"> <li>– English</li> <li>– Spanish</li> <li>– French</li> <li>– German</li> <li>– Italian</li> <li>– Russian</li> <li>– Japanese</li> <li>– Chinese</li> </ul> </li> <li>3. Check the communication settings. <ul style="list-style-type: none"> <li>– Port number is automatically detected after connection.</li> <li>– The PC software and balance have the same default communication Setting (see figure 3.2.2-1).</li> </ul> </li> <li>4. Click <b>OK</b> to continue.</li> </ol>
 <p style="text-align: center;">Figure 3.2.2-2</p>	<ol style="list-style-type: none"> <li>5. Once the balance has successfully connected to the PC, it shows the balance <b>Model</b> and <b>Serial Number</b> (see figure 3.2.2-2).</li> </ol>

## Manual Connection

	<ol style="list-style-type: none"> <li>1. Double-click the software.</li> <li>2. Choose the operation Language in the upper right corner.</li> <li>3. Click <b>Cancel</b> to continue.</li> </ol>
	<ol style="list-style-type: none"> <li>4. Choose the model from drop list, and then click <b>Enter</b>.</li> </ol>
	<ol style="list-style-type: none"> <li>5. Click <b>Setting - Connection</b> in the top.</li> <li>6. Check the connection setting as mentioned in <b>Automatic Connection</b>, and then click <b>OK</b> to confirm.</li> </ol>



Figure 3.2.2-3

7. Once the balance has successfully connected to the PC, it shows the balance **Model** and **Serial Number** (see figure 3.2.2-3).

### 3.2.3 Step 3: Customize Template

#### Main screen

After startup, OHAUS Label Designer appears with a blank state in the editing zone. The main screen has 4 sections.

- ❶ Import label template from PC / Export label template to PC in \*.tmp format
- ❷ Available values (print contents) read from balance
- ❸ Zone for editing label template or showing imported template
- ❹ Read label template from balance / Write label template to balance

Main screen when connecting to Explorer

## Customize label format

All label templates should be edited in ZPL programming language.

**Explorer: Maximum 100 lines with up to 80 characters in each line can be used in each label template.**

**Pioneer PX: Maximum 16 lines with up to 80 characters in each line can be used in each label template.**

To create a new label template, follow these steps:

	<p>1. Click <b>Read</b>, select Print Template, for example, #1 – 100mm * 75mm. Then click <b>OK</b> to continue.</p>
	<p>2. The below items can be customized by ZPL scripts</p> <ul style="list-style-type: none"> <li>– text position</li> <li>– font size</li> <li>– customized contents</li> <li>– Values read from balance (in blue text)</li> </ul>


**Example 1:**

To add “**date**” and “**time**” into the label, click **[5200] Date** and **[5201] Time** from the left side; and on the right side, ZPL scripts will show the code in blue text `[$5200!] [$5201!]` at the line.

Printout	ZPL Scripts		
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 5px;">Date/Time:</td> <td style="padding: 5px; text-align: center;">12/12/2021 05:39:23</td> </tr> </table>	Date/Time:	12/12/2021 05:39:23	<pre>^FT130,68^A0N,28,28^FH^FDDate/Time:^FS ^FT447,68^A0N,28,28^FH^FD[\$5200!] [\$5201!]^FS</pre>
Date/Time:	12/12/2021 05:39:23		

You can also define the font size of “date and time” and its position in the label.

- **28, 28** indicate the font size (height and width) of both “Date/Time” and `[$5200!] [$5201!]` (values read from Explorer).
- **Date/Time:** can be customized by users. For example, you can change it to **Date and Time:**.
- **130,68** indicate the position of **Date/Time:** in the label, i.e, from the left to right 130 dots, from top to down 68 dots. According to the label size, you can adjust the position, i.e, **150, 88**.

ZPL Scripts	^FT130,68^A0N,28,28^FH^FDDate/Time:^FS	^FT150,88^A0N,28,28^FH^FDDate and Time:^FS
<p><b>Printout</b></p>		



**Example 2:**

To print "**Sample ID**" as barcode, use below ZPL scripts and click **[5311] Sample ID** from left side, and on the right side, ZPL scripts will show the code at the line of "**^BY2,3,36^FT,447,129^BCN,,Y,N ^FD>\${5311}!^FS**".

Printout	ZPL Scripts
	<code>^BY2,3,36^FT,447,129^BCN,,Y,N ^FD&gt;\${5311}!^FS</code>

You can define the text position of bar code in the label.

- 447,129 indicate the position of **\${5311}!** (barcode) in the label, i.e, from the left to right 447 dots, from top to down 129 dots. According to the label size, you can adjust the position, i.e, 130, 129.

ZPL Scripts	Printout
<code>^BY2,3,36^FT,447,129^BCN,,Y,N ^FD&gt;\${5311}!^FS</code>	
<code>^BY2,3,36^FT,130,129^BCN,,Y,N ^FD&gt;\${5311}!^FS</code>	

To print "**Sample ID**" as number, use below ZPL scripts and click **[5311] Sample ID** from the left side, and on the right side, ZPL scripts will show the code at the line of "**^FT447,129^A0N,28,28^FH^FD\${5311}!^FS**".

You can define the font size of sample ID and the text position (see **Example 1**).

Printout	ZPL Scripts
<p>123456789012</p>	<code>^FT447,129^A0N,28,28^FH^FD\${5311}!^FS</code>

**Reset**

You can reset the edited label template to default template (#1 - 100mm\*75mm) by clicking **Reset**.

Model:



Serial Number:

Print Template:

### 3.2.4 Step 4: Write label template to balance


5 label templates (#1 - #5) are allowed to store in Explorer balances. Only one label template is allowed to store in Pioneer PX balances.

To write templates to the balance, follow below steps:

	<ul style="list-style-type: none"> <li>- Click <b>Write</b>.</li> <li>- Choose the template # you want to save to and define the template name. Click <b>OK</b> to continue.</li> </ul> <p>For example, if you choose #3, the original format in #3 will be replaced by new format automatically.</p>
	<ul style="list-style-type: none"> <li>- The left screen will appear if the template is written to balance successfully.</li> </ul>

To backup more templates to your PC, click **Export** and save the templates in format (\*.tmp).

These templates can be imported to OHAUS Label Designer for editing at any time.

For additional support for Zebra printers and label formats, contact Zebra directly at [www.zebra.com](http://www.zebra.com).

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 <p>The screenshot shows the OHAUS Label Designer interface. On the left, there is a tree view with categories like '重量' (Weight) and '模式' (Mode). The main area displays a list of templates with their corresponding codes and descriptions. A '选择模板' (Select Template) dialog box is open, showing '编号' (Number) set to '#3' and '名称' (Name) set to '3'. At the bottom, there are fields for '型号' (Model: EX225D) and '序列号' (Serial Number: T1950930885), along with buttons for '打印模板' (Print Template), '读取' (Read), and '写入' (Write).</p>	<ul style="list-style-type: none"> <li>- 点击<b>写入</b>。</li> <li>- 选择所需写入的模板编号# 并定义模板名称。点击 <b>OK</b> 继续。</li> </ul> <p>例如，如果您选择模板#3,那么模板#3 中的原始标签格式将自动替换为新存入的标签格式。</p>
 <p>The screenshot shows an '信息' (Information) dialog box with a blue information icon and the text '写入数据到天平成功。' (Data written to the scale successfully.) Below the message is an 'OK' button.</p>	<ul style="list-style-type: none"> <li>- 如果模板成功写入天平，则会显示左图。</li> </ul>
 <p>The screenshot shows the 'RS232接口设置' (RS232 Interface Settings) menu. It includes options for '通讯传输设置' (Communication Transfer Settings), '打印设置' (Print Settings), '标签打印' (Label Print) which is '开启' (On), and '标签模板' (Label Template) which is set to '3'. At the bottom, there are buttons for '主界面' (Main Interface), '后退' (Back), and '完成' (Finish).</p>	<ul style="list-style-type: none"> <li>- 连接天平和打印机，选择模板“3”后回到称量界面，按打印键，打印。</li> </ul> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;"> <p><b>生产过程称量记录</b></p> <p>物料名称: ABC</p> <p>物料批号: 12345678</p> <p>生产批号: C234567</p> <p>毛重: 10.000g</p> <p>净重: 10.000g</p> <p>皮重: 0.000g</p> <p>称量人: User 1</p> <p>复核人:</p> <p>称量日期: 12/12/2021 05:39:23</p> <p>备注:</p> </div>

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P/N 30708369C © 2022 Ohaus Corporation, all rights reserved / todos los derechos reservados / tous droits réservés
Ohaus Corporation