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Ziegelei 1
D-72336 Balingen
E-Mail: info@kern-sohn.com

Phone: +49-[0]7433- 9933-0
Fax: +49-[0]7433-9933-149
Internet: www.kern-sohn.com

Operating instructions

KERN App Count - Piece-count function

KERN SET-032

Version 1.0
2020-06
GB

The current version of these instructions can also be found online under:
<https://www.kern-sohn.com/shop/de/DOWNLOADS/>
Under the section Operating instructions

SET-032-BA-e-2010_count

KERN App Count
Version 1.0 2020-06
Operating instructions
SET-032


Contents

1	Weighing function Count – Piece-count function	3
1.1	Determining the reference weight	3
1.2	Determined weight shall correspond to 5 items.....	5
1.3	Certain weight has to correspond to a certain quantity of items.....	7
1.4	Save the result data	8

1 Weighing function Count – Piece-count function

The piece-count function helps to determine the piece quantity of a number of items, referring to a previously defined reference weight.

In the menu FUNCTIONS click button „Quantity“. The start screen for piece counting appears.

1.1 Determining the reference weight

In order to determine the reference weight, click the button „Define reference weight“. The screen to determine the reference weight appears. The reference weight can either be determined by weighing or by entering the weight value.

By weighing:

Place the sample which corresponds to the piece quantity „1,“ on the weighing plate and click the button „Apply weight“:

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Click CONFIRM. The piece quantity „1,, is now displayed, as well as the reference weight and the total weight:

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
Now you can weigh the parts. Remove the sample and place the parts to be counted, the quantity, as well as the reference weight and the total weight are displayed:

1.2 Determined weight shall correspond to 5 items

Place a weight which shall correspond to 5 items, on the weighing plate and press



button :

The piece quantity is displayed. The reference weight has also been calculated automatically and the total weight is displayed.

To count unload the weighing plate and place the parts upon. The quantity of the placed parts is displayed, as well as the reference weight and the total weight:

The weight for 10 items  and 20 items  can also be determined.

1.3 Certain weight has to correspond to a certain quantity of items

If you want to allocate a self-defined weight to a self-defined quantity of items, place

a sample of that weight and press button . Now you can enter the quantity of items which shall match with this weight:

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Referenzgewicht *
123,95

Einheit
(g) ▾

Then click CONFIRM, reference weight and total weight are displayed.

Remove the sample and place the quantity of parts to be counted:

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TARA 0,0 kg 0

Referenzgewicht 123,90 g / St. 📄

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Same procedure with all the other parts.

1.4 Save the result data

To save the weighing result, press button

Enter ID and name. The result can be saved as master data object and printed out.

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Fa. Kern und Sohn GmbH